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EN Instructions for use

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NO Brukerveiledning

SV Bruksanvisning

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SK Návod na používanie

HR Upute za uporabu

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1. Device identification

MD The SQY Foot is a class I medical device

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Phone: +33 (0)2 28 23 67 33

Email address: contact@sqyinnovations.com


Website: www.sqyinnovations.com

2. Device description

The SQY Foot is a dynamic energy-return prosthetic foot composed of:

- A. A split lower carbon fiber blade
- B. An upper carbon fiber blade
- C. A titanium pyramid adapter
- D. Drainage grooves
- E. Titanium screws
- F. An encapsulating cosmetic shell made of elastomer

The device is hereinafter referred to as “the device.”

3. Intended use and indications

This device is intended to be integrated into a custom-made external lower limb prosthesis to provide foot function for patients with unilateral or bilateral amputation and/or congenital lower limb deficiency. The device is indicated for patients with a moderate to high activity level (K3), for walking and physical activities without excessive overload.

Maximum user weight (including load): 146 kg.

The device is supplied exclusively to healthcare professionals (prosthetists). The prescription is established by a physician. The device is intended for multiple use on a single patient.

4. Properties

The device has been tested in accordance with **ISO 10328** for a maximum user weight of 146 kg, for 2 million cycles.

Service life depends on activity level, conditions of use, and compliance with maintenance instructions.

Selection of foot category according to patient weight and impact level

Weight (kg)	Impact level		
	Low	Medium	High
44-51	1	1	2
52-58	1	2	3
59-67	2	3	4
68-76	3	4	5
77-87	4	5	6
88-99	5	6	7
100-115	6	7	8
116-129	7	8	9
130-146	8	9	/

Low: Daily walking and occasional sports such as golf, cycling

Medium: Sustained walking, frequent or daily sports such as jogging, tennis, volleyball

High: Daily activities such as running, climbing, lifting and carrying heavy objects, football.

5. Mechanism of Action

The long lower sole of the foot allows for a continuous and smooth rollover during walking. The lower sole of the foot is made of a split carbon blade running lengthwise. The carbon provides the foot with great flexibility and lightness. The upper carbon blade gives the foot the ability to store and return energy. The combination of the two split blades allows for inversion/eversion of the foot.

The cosmetic cover of the foot encapsulates the carbon fiber blade and provides additional protection as well as shock absorption for the user.

6. Intended Use / Indications

This medical device is supplied to healthcare professionals (prosthetists) responsible for training the patient in its use.

The prescription is established by a physician who assesses the patient's ability to use the device. This device is intended for multiple use on a single patient. It must not be used on another patient.

This device is intended to be integrated into a custom-made external lower limb prosthesis to provide foot function for patients with unilateral or bilateral amputation and/or congenital lower limb deficiency. This device is indicated for patients with a moderate to high activity level (K2 to K4) for walking and physical activities without excessive overload.

Maximum weight (including load carried): **146 kg**.

Warning: Do not exceed the weight limit. Risk of prosthetic foot failure. Incorrect category selection may also result in improper functioning of the device. This device is not recommended for sports competitions.

7. Clinical Benefits

The device ensures optimal comfort and safety for transtibial or transfemoral amputees during movement, on uneven terrain, slopes, and stairs. The SQY Foot also provides the following benefits:

- Walking comfort
- Adaptation to uneven terrain
- Slip resistance on wet surfaces
- Shock absorption and noise reduction

8. Accessories and Compatibility

The foot includes a male pyramid attachment designed to be compatible with standard female pyramid connectors.

The SQY Foot comes with an attachment accessory for open shoes such as sandals and flip-flops.

9. Alignment

Basic Alignment

Note: Prioritize knee alignment over foot alignment in case of discrepancy.

Sagittal Plane


Align the socket in flexion and the heel in height according to the needs of the amputee.

Position the load line, from the center of the socket at the patellar tendon, so that it corresponds to the anterior edge of the pylon.

Frontal Plane

Align the socket in adduction/abduction according to the needs of the amputee.

Position the load line, from the center of the socket at the patellar tendon, so that it passes through the midline of the foot module in the neutral M-L position.

Static Alignment

- Ensure that the patient stands with weight evenly distributed on both legs.
- Check the appropriate length of the prosthesis.
- Check internal/external rotation.
- Check that the load on the toe and heel is correct.

Dynamic Alignment

The rollover and propulsion can be influenced by the following factors:

- The anterior/posterior position of the device.
- Plantar/dorsiflexion.
- Shoe characteristics.

Consider the following actions if necessary:

Symptoms	Actions
<ul style="list-style-type: none"> • The device comes to flat too early (the patient feels like they are falling into a hole) • Additional effort is required to roll over the toe • The toe is perceived as too stiff • The knee goes into hyperextension 	<ul style="list-style-type: none"> • Perform an anterior translation of the socket (or a posterior translation of the device) • Consider dorsiflexion • Check the shoe heel and shoe performance
<ul style="list-style-type: none"> • The movement from heel to toe is too fast. • Poor prosthesis control at initial contact • Very little energy return is felt • Toe-off from the toe is too weak • The knee becomes unstable 	<ul style="list-style-type: none"> • Perform an anterior translation of the socket (or a posterior translation of the device) • Consider plantarflexion • Check the shoe heel and shoe performance

10. Detection of malfunctions

Any serious incident involving the device must be reported to the manufacturer and the competent authorities.

11. Warning

Never loosen the pyramid fixation screws.

The patient must inform their prosthetist if they gain or lose a significant amount of weight.

If the patient notices abnormal behavior or experiences changes in the characteristics of the device (noise, play, excessive wear, etc.), or if the device has suffered a significant impact, they must stop using the device and consult their prosthetist.

Failure to follow the instructions for use is dangerous and will void the warranty.

The use of components from other manufacturers may result in overloading the device.

The device has a positive buoyancy.

Use appropriate shoes that are compatible with the foot cosmesis.

Any serious incident must be reported to the manufacturer and the competent authority.

The prosthetic foot must always be used with suitable shoes on abrasive or sharp surfaces.

12. Contraindications

No known contraindications when the device is used in accordance with the specified indications and limitations.

13. Side Effects

No side effects have been reported to date.

Any serious incident occurring in connection with the device must be reported to the manufacturer and the competent authority of the Member State in which the user is established.

14. Maintenance and Inspection

The device and the prosthesis as a whole must be examined by a healthcare professional. The interval should be determined according to the patient's activity level.

No lubrication, screw adjustment, or other intervention is required for this foot except by the prosthetist.

15. Cleaning

After use in water, sand, or a dusty environment:

- Rinse the foot with clean water
- Dry thoroughly with a soft cloth

DO NOT use aggressive detergents or solvents.

16. Environmental conditions

The device can tolerate contact with: salt water, chlorinated water, and mild soaps. It can also tolerate occasional exposure to sand, dust, and dirt. Continuous exposure is not permitted. Dry with a cloth after contact with fresh water or moisture. Clean with fresh water after accidental exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

17. Disposal

The device and packaging must be disposed of in accordance with applicable local or national environmental regulations.

18. Regulatory Information

This product is a CE-marked medical device and certified in compliance with Regulation (EU) 2017/745.

Device identification information is provided on the packaging.

19. General Safety Instructions

WARNING: The use of a lower limb prosthesis involves an inherent risk of falling, which may result in injury. The healthcare professional must provide the patient with all the information contained in this document necessary for the safe use of this device.

WARNING: If the functionality of the device changes or decreases, or if the device shows signs of damage or wear preventing it from functioning normally, the patient must stop using it and contact a healthcare professional.

WARNING: Risk of structural failure. Components from other manufacturers have not been tested and may result in excessive load on the device.

WARNING: Ensure proper fixation by following the applicable device assembly instructions.

WARNING: Improper use of the foot, contrary to the recommendations of your prosthetist, may cause degradation of foot components (for example, carrying heavy loads, excessive stress, exceeding the service life, etc.).

WARNING: Always use a ramp when going down stairs and at any other time if possible.

WARNING: Appropriate shoes, well fitted to the foot cosmesis, must be used.

WARNING: Ground grip – Unlike conventional prosthetic foot cosmesis, the SQY Foot is made from a material with increased ground grip. This feature may alter the usual sensations during the first uses. A gradual adaptation phase is recommended.

WARNING: Use of the prosthetic foot on rough surfaces (stones, rocks) may damage the foot and reduce its grip.

WARNING: Always use appropriate footwear on rough or sharp surfaces.

WARNING: This foot has a positive buoyancy and this must be taken into account when calculating the overall buoyancy of the prosthesis.

20. Liability

The reliability and compatibility of the SQY Foot prosthetic device with custom prosthetic sockets, when used as intended, are verified. AQUALEG SAS declines all liability in the following cases:

The device is not maintained as indicated in the instructions for use.

The device is not assembled according to the manufacturers' recommendations.

The device is used outside the recommended scope regarding conditions of use, application, or environment.

21. Description of the symbols

	Medical Device		Manufacturer		CE Marking
	Warning		Single patient, multiple use		eIFU for patients

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